



*Arglaes provides a seven-day, non-cytotoxic barrier against infection*



**Arglaes® Controlled-Release Silver Technology**

# Reduce bioburden with Arglaes® Antimicrobial Silver Barrier Dressings



Arglaes began the antimicrobial silver revolution as the first product to provide controlled-release, ionic silver for up to seven days.

Arglaes technology utilizes ionic silver to create an environment hostile to bacteria and fungi, but completely non-cytotoxic. The sustained-activity ionic silver maintains full efficacy for up to seven days. The Arglaes technology is available as a transparent film dressing, with and without an alginate pad. The film has a high moisture vapor transmission rate. Arglaes Powder is well suited for difficult-to-dress wounds as its alginate base forms a soft, conforming gel when it contacts wound exudate.

## **Reduces bioburden**

Effective against a broad spectrum of bacteria and fungi (including MRSA, VRE and *E.coli*).

## **Constant antimicrobial protection**

Controlled-release antimicrobial silver.

## **Non-cytotoxic**

Controlled-release polymers deliver small, non-cytotoxic amounts of silver that kill bacteria and fungi, but do not harm healthy tissue.

## **Extended wear time**

Effective for up to seven days.

## **Reduce costs**

Dramatically reduces costs associated with treating infections.

## **No known resistance**

Ionic silver has shown no known resistance to any bacteria or fungi, including antibiotic resistant bacteria.



## The Arglaes family

---



*Arglaes Film*

### **Ideal for post-op and line sites**

Arglaes Film is perfect for managing bioburden on line sites, post-operative incisions and donor sites. Up to seven-day wear time



*Arglaes Island*

### **Manages fluid and bioburden**

Arglaes Island features a calcium alginate pad for fluid management in addition to controlled-release silver. Up to 5-day wear time.



*Arglaes Powder*


### **Ideal for deep, tunneling and highly exuding wounds**

Arglaes Powder combines controlled-release silver with an alginate powder to offer bioburden reduction and fluid management to virtually any size, shape or depth of wound. Up to five-day wear time.

## Arglaes in wound care

---

*Arglaes Powder delivers controlled-release antimicrobial silver to any size, shape or depth of wound.*



As Arglaes Powder mixes with wound exudate, it turns into a gel that adheres to the wound bed

Any size. Any shape. Any depth. No matter what the wound<sup>1</sup>, use Arglaes Powder and you'll reduce the risk of infection.

Utilizing controlled-release polymers, Arglaes Powder delivers a constant stream of antimicrobial silver ions into the wound over a period of five days.

Continuous delivery, at a constant rate, means that only minute quantities of silver ions are required to maintain a continuous antimicrobial barrier without cytotoxicity.

In addition, Arglaes Powder contains alginate to aid in fluid handling. As the powder mixes with wound exudate, it turns into a gel that adheres to the wound bed and is easily removed during wound irrigation.

Arglaes Powder is easily combined with other dressings to create a system for bioburden control and optimal moist wound healing.

<sup>1</sup> See indications

## Arglaes in the operating room

---

### *Arglaes Film and Arglaes Island offer post-operative infection control.*

Post-operative infections are obviously a major problem in terms of trauma to the patient and the costs associated with treating the infection. Arglaes can help protect incision sites by providing up to a seven-day barrier against infection.

#### **Arglaes is indicated for a variety of post-operative settings including:**

- Sternotomy
- CABG procedures
- FEMPOP procedures
- Open-heart procedures
- Long-term catheters
- Skin grafts
- Donor sites
- Tram flap

“We had about a 1½% rate of sternal wound infections. Since we started using the Arglaes dressings, we have not had a deep sternal infection. Arglaes has significantly reduced our complication rate and I would recommend it for any program doing cardiac surgery.”

Edward Pottmeyer, M.D.  
Cardiothoracic Surgeon , Mercy Medical Center, Redding, CA



Arglaes Island applied on a sternal incision.



Arglaes Film (4" x 4¾") applied on a minimally invasive incision for saphenous vein harvesting.




Arglaes Island (4" x 4¾") applied on a minimally invasive sternal incision.

## Arglaes on line sites

---

### *Arglaes Film helps provide a barrier against infection on line sites.*

Infections on line sites are a costly problem to treat. Arglaes Film provides an easy-to-use and effective barrier against line site infections. Simply apply as you would a transparent film, and Arglaes delivers controlled-release silver, inhibiting the growth of bacteria and fungi, for up to seven days.



Arglaes Film applied on a line site.

# Arglaes Clinical Information

## Zones of Inhibition from Sustained Antimicrobial Activity

Experiments were conducted to assess the antimicrobial activity of a controlled-release silver compound against a range of microbiological species. Each species was selected because of its significance as a potential pathogen.

Agar plates inoculated with the test organisms were prepared and 10mm discs of the controlled-release silver compound and a known positive control of povidone iodine were placed on a series of separate plates.

Zones of inhibition were recorded at 24 hours, 48 hours and 7 days.

## Arglaes Zone of Inhibition Study (in mm)

Organism	24 hours			48 hours			7 days		
	Arglaes	Iodine (Povidone Iodine)	Fertility Control	Arglaes	Iodine (Povidone Iodine)	Fertility Control	Arglaes	Iodine (Povidone Iodine)	Fertility Control
<i>Proteus vulgaris</i>	11.3	12.1	+++	10.2	5.2	+++	9.3	2.9	+++
<i>Acinetobacter baumannii</i>	14.8	6.8	+++	10.2	4.0	+++	10.5	3.8	+++
<i>Enterococci faecium</i>	9.2	0.0	+++	8.8	0.0	+++	9.5	1.2	+++
<i>Serratia marcescens</i>	9.0	4.4	+++	7.5	0.0	+++	7.3	0.0	+++
<i>Candida albicans</i>	13.5	15.1	+++	13.1	8.2	+++	15.6	3.8	+++
<i>Pseudomonas aeruginosa</i>	9.7	1.4	+++	8.9	0.0	+++	8.0	0.0	+++
<i>Proteus mirabilis</i>	8.8	6.9	+++	8.1	5.7	+++	9.9	5.9	+++
<i>Staphylococcus aureus</i>	14.7	0.0	+++	13.5	0.0	+++	13.4	0.0	+++
<i>Escherichia coli</i>	7.5	0.7	+++	6.9	0.6	+++	6.7	0.0	+++
<i>Enterobacter cloacae</i>	5.4	8.0	+++	5.1	0.0	+++	4.8	0.0	+++
<i>Staphylococcus aureus</i> (MRSA, NCTC 12493*)	20.0	10.6	+++	20.5	0.0	+++	20.3	0.0	+++
<i>Klebsiella Edwardsii</i> var <i>Edwardsii</i>	11.3	12.1	+++	10.2	5.2	+++	9.3	2.9	+++
<i>Staphylococcus aureus</i> (MRSA, NCTC 12232*)	13.9	2.6	+++	12.8	1.7	+++	12.5	1.5	+++
<i>Staphylococcus epidermidis</i>	16.0	0.0	+++	15.7	0.0	+++	22.1	0.0	+++

\*National Collection Type Culture Catalog. Independent study performed by Wickham Laboratories Limited, Hampshire, England




## Antimicrobial Activity

Six separate experiments have been completed to test the in vitro activity of controlled-release barrier film dressing (Arglaes) and adhesive, which also contains controlled-release compound.

Results from all six studies confirmed the antimicrobial activity of the controlled-release barrier film dressing against the pathogens commonly associated with wound infections.

Pathogens tested include: *Staphylococcus aureus*, *Escherichia coli*, *Faecal streptococcus*.




## Sustained-Release Activity of Arglaes



Independent study performed by Wickham Laboratories Limited, Hampshire, England

## Illustration of Controlled-Release Polymers

Arglaes provides an antimicrobial barrier for up to seven days. Silver ions are released at a constant rate, providing a barrier that inhibits the growth of new bacteria and helps prevent the migration of already existing bacteria.



Illustrations shown here are representations of sustained-activity antimicrobial barrier effect.

# Arglaes Product Information

## Indications

- Pressure ulcers (Stage I-IV)
- Partial and full thickness wounds
- Leg ulcers
- Diabetic ulcers
- Central lines, CVPs and PICC lines (Arglaes Film only)
- Surgical wounds
- Grafted wounds (Arglaes Powder only)
- Donor sites
- Lacerations and abrasions
- 1st and 2nd degree burns

## Contraindications

- Third degree burns
- Individuals with a known sensitivity to silver
- As a surgical implant

## Change Frequency

- Arglaes Film may be left in place for up to 7 days.
- Arglaes Island and Arglaes Powder may be left in place for up to 5 days.
- Dressing change frequency will depend upon the amount of exudate.



## Arglaes Ordering Information

<u>Item No.</u>	<u>Description</u>	<u>Packaging</u>
MSC9023	Film, 2 $\frac{3}{8}$ " x 3 $\frac{1}{8}$ "	10/bx, 100/cs
MSC9045	Film, 4" x 4 $\frac{3}{4}$ "	10/bx, 100/cs
MSC9069	Film, 4 $\frac{3}{4}$ " x 10"	10/bx, 50/cs
MSC9314	Film, 3" x 14", Post-Op Style	10/bx, 50/cs
MSC9123	Alginate Island, 2 $\frac{3}{8}$ " x 3 $\frac{1}{8}$ ", 1" x 2" pad	10/bx, 100/cs
MSC9145	Alginate Island, 4" x 4 $\frac{3}{4}$ ", 2" x 2" pad	10/bx, 100/cs
MSC9169	Alginate Island, 4 $\frac{3}{4}$ " x 10", 2 $\frac{3}{4}$ " x 8" pad	10/bx, 50/cs
MSC9210	Powder, 10 gm bottle	5/bx, 20/cs
MSC9210SP	Powder, 10 gm bottle in Sterile Pouch	5/bx, 20/cs
MSC9205	Power, 5 gm bottle	5/bx, 20/cs

1-800-MEDLINE | [www.medline.com](http://www.medline.com)

© 2003 Medline Industries, Inc., Mundelein, Illinois 60060

Medline is a registered trademark of Medline Industries, Inc. Arglaes is a registered trademark of Giltech, Ltd.

MKT2033/LIT034/20M