



# Medline Urologicals

## Better By Design



# Medline Urologicals...Better By Design

## Silvertouch™ Foley Catheters

Latex-free catheters incorporate the power of silver through a patented process that binds silver ions to the catheter's lubricious coating. The hydrophilic coating is covalently bonded to the catheter surface, hydrates quickly, and maintains its lubricity for at least seven days, helping to reduce friction and irritation.



## Pre-Saturated Swabsticks

Swabsticks offer uniform, complete saturation and encourage good technique by keeping hands away from the patient. Swabsticks are wound tightly and pre-saturation assures rayon fibers will not fly through air and settle on patient or catheter.



## Needleless Sample Port

Medline was the first to offer a needleless sample port, providing a safe way to take urine samples. Made from self sealing thermo-plastic elastomer that is designed to be used with a blunt cannula or sharp needle.



## Ultimate Anti-Reflux Device

Delivers the highest quality reflux control with a mylar valve which guarantees one-way flow. Mylar valve is easily seen through the clear drip chamber, allowing visualization of urine flow. Ideal for the hospital environment.



## Anti-Reflux Tower

Specially designed to provide reflux protection without mechanical parts to hinder flow. Features a non-wettable filter which helps prevent vacuum action and facilitates drainage. Designed to meet the needs of the extended care environment.



## Slide-Tap

Available from Medline, the Slide-Tap Mechanism allows for easier one hand bag drainage.



## Metal Clamp

Designed to accommodate all health care practitioners, as well as patients with limited hand dexterity.



## Double Hanger

Redesigned double hanger made of reinforced flexible PVC can easily rest on patient's bed or wheelchair. Ribbed frame provides extra support for full or heavy bags.



## STAR™ Tubing

Exclusive to Medline, STAR Tubing is the first and only urinary drain bag tubing equipped with ribbed inner lumen to enhance urine outflow, including thick or fibrous urine. Reinforced tubing helps to prevent kinking.



## Aloetouch® 3G Gloves

Medline's innovative latex-free, powder-free stretch synthetic exam gloves give outstanding tactile sensitivity with a unique softness and flexibility, and coated with Aloe Vera to help moisturize and soothe dry hands.



## Table of Contents

Catheter Associated Urinary Tract Infections	4-5	Foley Catheter Accessories	16
Proper Drainage	6	External Catheters	17
CDC Guidelines	7	Infant/Pediatric Products	17
Silvertouch Foley Catheters	8-9	Intermittent Catheters & Trays	18-19
Closed System Foley Trays	10-11	Specimen Collection & Containers	20-22
Foley Catheter Insertion Trays & Leg Bags	12	Irrigation Syringes & Trays	23
Urinary Drainage Bags & Urine Meters	13	Product Cross Reference	24-25
Foley Catheters	14-15	Index	26-27

# Catheter Associated Urinary Tract Infections

## WHO'S AT RISK?

Any patient with a urinary catheter is at risk for a urinary tract infection (UTI). One of the most important risk factors for UTI is duration of catheterization. Within hours, a biofilm develops on the surface of an indwelling catheter and the drainage system. This biofilm has been linked to persistent bacteriuria; as the duration of catheterization increases, so does the prevalence of bacteriuria and risk of infection.<sup>1</sup>

### Symptoms of UTI

Every bedside nurse should recognize these common symptoms of UTI:

- Back or lower abdominal pain
- Fever
- Cloudy or bloody urine
- Mucous discharge
- Suprapubic pain
- Worsening mental function

### Treatment for Symptomatic Patients

Many catheterized patients will develop bacteriuria but never display symptoms of UTI. Antibiotic treatment is recommended only in the presence of symptoms. A catheter should be replaced in a patient presenting symptoms, and the sample drawn from the new catheter. Removing the biofilm-laden catheter from the patient assures that any pathogens cultured came from the urinary tract, not the previous catheter.<sup>2</sup>

### Pathogens Associated with UTI

- Escherichia coli
- Klebsiella
- Candida albicans
- Enterococcus
- Pseudomonas
- Proteus
- Enterobacter
- Staphylococcus

### CAUTI Facts:

- The most common nosocomial infection
- The second most common cause of nosocomial bloodstream infection
- UTIs account for 40% of all nosocomial infections
- 25% of hospitalized patients receive an indwelling catheter
- Catheterized patients have a 5% daily risk of bacteriuria
- 10–25% of patients with bacteriuria will have symptomatic UTI

#### References:

- 1 Saint S, Chenoweth CE. Biofilms and catheter-associated urinary tract infections. *Infectious Disease Clinics in North America*, Jun 2003;17(2):411-32.
- 2 Raz R, Schiller D, Nicolle LE. Replacement of catheter improves the outcome of patients with permanent urinary catheter and symptomatic bacteriuria [abstract]. Presented at the 38th Annual Interscience Conference on Antimicrobial Agents and Chemotherapy (ICAAC); 1998 Sep 24–27; San Diego, CA: 532.

#### Other sources:

Guideline for Prevention of Catheter-associated Urinary Tract Infections. Centers for Disease Control and Prevention. [www.cdc.gov](http://www.cdc.gov).  
 Society of Urological Nurses and Associates, [www.suna.org](http://www.suna.org).  
 Guideline for Hand Hygiene in Health-Care Settings. Centers for Disease Control and Prevention, 2002. [www.cdc.gov](http://www.cdc.gov).

## Education—the BEST Prevention Plan

Education is a key factor in CAUTI prevention. Make sure that anyone catheterizing a patient is educated on the proper technique for catheter insertion and care. Review of proper

technique should occur periodically and should include the concepts discussed in this article as well as any updated clinical guidelines.

## EIGHT STEPS TO PREVENTION:

- 1 Handwashing.** This remains the single most important step in preventing the spread of infection. Hands should be cleansed immediately before and after any contact with or manipulation of the catheter site, tubing, or drainage bag.
- 2 Catheterize only when necessary.** Acute urinary retention or obstruction and the need to observe urinary output are the most common reasons for Foley catheterization in hospitals. Most in-house catheters are initially placed in the operating room or emergency department. Review your facility's policies to determine if they clearly state indications for catheter placement and call for removal as soon as clinically appropriate. Other methods of urinary drainage (intermittent catheters or male external catheters) should be used whenever possible. Be sure to maintain appropriate documentation and communication when patients are transferred.
- 3 Aseptic preparation and insertion.** In addition to diligent handwashing, sterile technique should be followed when inserting the catheter. Check for Foley balloon patency and apply lubricating jelly to the tip of the catheter prior to prepping the skin. Once a sterile-gloved hand touches the patient's skin, it should be considered contaminated. Any sterile equipment that contacts the patient should not re-enter the sterile field.
- 4 Closed system.** The catheter should remain connected to the closed drainage system. Breaking the connection can create an entry point for bacteria. Empty the drainage bag regularly, taking care not to touch the spigot to the nonsterile collecting container.
- 5 Good drainage.** Remember “Bag Below the Bladder” so gravity drains urine away from the bladder. Keep the tubing as straight as possible, and avoid kinking or compression. To avoid urine reflux and saturation of the vent, a drainage bag should not be left horizontal on any flat surface, especially face-down.
- 6 Smallest catheter possible.** To minimize urethral trauma, use the smallest catheter that also allows for good drainage.
- 7 Secure the catheter.** To prevent movement of the catheter and decrease the risk of bleeding, trauma, necrosis, and bladder spasms from pressure and traction, the catheter should be secured to the patient's thigh or abdomen after insertion.
- 8 Perineal cleansing.** A regular meatal cleansing regimen is not recommended for UTI prevention. However, daily perineal skin care is still important. Take extra care in cleaning patients after every incontinent stool episode.

## Silver: Another Weapon Against Bacteria

Because silver is a broad-spectrum antimicrobial not associated with bacterial resistance, catheters coated with ionic silver hydrogel have gained considerable attention and acceptance. Silver is effective against gram-negative and gram-positive bacteria as well as fungi. Even antibiotic-resistant bacteria (such as MRSA and VRE) are vulnerable to ionic silver. The combination of silver and a high-quality hydrogel can facilitate the controlled release of silver (on contact with body fluids) over several days.

## Closed System Drainage

A closed urinary drainage system reduces the risk of catheter-associated urinary tract infection; any disconnection in the closed system allows for potential introduction of microorganisms to the urinary tract.

*From CDC Guideline for Prevention of Catheter-associated Urinary Tract Infection:*

### Closed Sterile Drainage

- a. A sterile, continuously closed drainage system should be maintained. Category I
- b. The catheter and drainage tube should not be disconnected unless the catheter must be irrigated. Category I

## Aseptic Insertion Technique

A urinary catheter will pass through the urethra and into the bladder; in order to minimize risk of infection associated with insertion, aseptic technique should always be maintained during the insertion process.

*From CDC Guideline for Prevention of Catheter-associated Urinary Tract Infection:*

### Catheter Insertion

- a. Catheters should be inserted using aseptic technique and sterile equipment. Category I
- b. Gloves, drape, sponges, an appropriate antiseptic solution for periurethral cleaning, and a single-use packet of lubricant jelly should be used for insertion. Category II

## Good Drainage

One primary infection control objective in urinary drainage is to maintain flow. By keeping urine flowing away from the patient's bladder, you can minimize the risk of bacteriuria in the patient's urinary tract that could lead to UTI.

### In order to maintain flow:

- Drainage bag should hang below the level of the patient's bladder.
- Keep the tubing as straight as possible from catheter to drainage bag.
- Avoid kinking or compression of the drainage tubing.
- Medline utilizes STAR Tubing (extruded with multiple ridges on the interior), which is very difficult to obstruct drainage by accidental compression. Also, while a kink can be manipulated into the tubing, the STAR design encourages a prompt spring-back, so the kink is not maintained.
- Avoid reflux from drainage bag into the tubing.
- Consider continuous irrigation when necessary.
- Secure the catheter. To minimize catheter movement and kinking, secure the catheter with a securement device or catheter strap, or alternatively, with tape.



**Closed System**  
Sealed to maintain closed system



**Aseptic Technique**  
Sterile equipment



**Good Drainage**  
STAR Tubing

## CDC Guidelines for Prevention of Catheter-Associated Urinary Tract Infections

### Summary of Major Recommendations

#### Category I. Strongly Recommended for Adoption

- Educate personnel in correct techniques of catheter insertion and care.
- Catheterize only when necessary.
- Emphasize handwashing.
- Insert catheter using aseptic technique and sterile equipment.
- Secure catheter properly.
- Maintain closed sterile drainage.
- Obtain urine samples aseptically.
- Maintain unobstructed urine flow.

#### Category II. Moderately Recommended for Adoption

- Periodically re-educate personnel in catheter care.
- Use smallest suitable bore catheter.
- Avoid irrigation unless needed to prevent or relieve obstruction.
- Refrain from daily meatal care with either of the regimens discussed in text.
- Do not change catheters at arbitrary fixed intervals.

#### Category III. Weakly Recommended for Adoption

- Consider alternative techniques of urinary drainage before using an indwelling urethral catheter.
- Replace the collecting system when sterile closed drainage has been violated.
- Spatially separate infected and uninfected patients with indwelling catheters.
- Avoid routine bacteriologic monitoring.



Silvertouch. The Standard for Silver.

**Silver ions provide antimicrobial action**

Silver ions have long been recognized for their broad-spectrum antimicrobial action. Ionic silver is effective against gram-positive and gram-negative bacteria, as well as fungi. Ionic silver hydrophilic catheter coatings can be very effective against the most common and serious UTI pathogens, including E. coli, Candida albicans, Staph aureus, Pseudomonas, Klebsiella, MRSA, VRE, and others.\* Silvertouch catheters incorporate the power of silver through a patented process that binds silver ions to the catheter's lubricious coating.

**Hydrophilic coating ensures a week's worth of comfort**

Every Silvertouch catheter has a coating that hydrates rapidly and maintains its lubricity for at least seven days, helping to reduce friction and irritation. The hydrophilic coating is covalently bonded to the silicone catheter surface and will not peel, flake, or crack.

**Look for the green!**

Why are Silvertouch catheters green? The brilliant color is visual evidence of the silver hydrophilic coating on each catheter.

**Latex-Free**

Silvertouch catheters are latex-free, offering patient and caregiver safety, with the advantages of a 100% silicone catheter.

\*In-Vitro test data on file.



Silvertouch Foley Tray



## Did You Know?

- Catheter-Associated UTIs (CAUTI) represent approximately 40% of all health care-acquired infections
- 25% of hospitalized patients receive an indwelling catheter
- 5% daily risk of CAUTI with short-term catheterization
- 10-25% of patients with bacteriuria will have symptomatic UTI
- A urinary tract infection adds 1-4 days to a hospital stay at an annual cost of \$500 million.

### Ordering Information

Item Number	Description	Pkg
<b>Silvertouch Foley Catheters</b>		
DYND141014	14Fr 5cc	10 ea/cs
DYND141016	16Fr 5cc	10 ea/cs
DYND141018	18Fr 5cc	10 ea/cs
DYND141020	20Fr 5cc	10 ea/cs
DYND141216	16Fr 30cc	10 ea/cs
DYND141218	18Fr 30cc	10 ea/cs
DYND141220	20Fr 30cc	10 ea/cs

### Silvertouch Foley Catheter Trays

DYND140214	14Fr 5cc	10 ea/cs
DYND140216	16Fr 5cc	10 ea/cs
DYND140218	18Fr 5cc	10 ea/cs

### Silvertouch Foley Urine Meter Trays

DYND140414	14Fr 5cc	10 ea/cs
DYND140416	16Fr 5cc	10 ea/cs
DYND140418	18Fr 5cc	10 ea/cs

510k Test Data available upon request.

**“The best way to reduce the morbidity, mortality, and associated costs of catheter-associated UTI is to prevent or minimize bacteriuria.”**

(Ref: Alcaid & Lichtstein, Hospital Physician, August 2004)

**“Catheters coated with silver alloy/hydrogel have recently been introduced into practice, and a growing body of literature supports their use in a variety of clinical settings.”**

(Ref: Rupp et al; AJIC December 2004)

## Closed System Foley Trays

- Sealed junction insures maximum protection against infection while securing Foley catheter to drainage tubing.
- Urine Bag Drainage Port has a Slide-Tap valve with "see through" drain tube.
- Needleless Sample Port allows urine sampling with standard luer slip syringe or blunt cannula.
- The Medline Prep Tray includes: Aloetouch 3G latex-free wrapped gloves, 3 PVP swabsticks, underpad, fenestrated drape, prefilled inflation syringe, lubricant syringe, specimen container w/ ID label.
- Exclusively from Medline – Aloetouch 3G powder-free stretch synthetic exam gloves.

### Medline



#### Clean Uniform Saturation

Provide better aseptic technique  
Tightly wound fibers

### Others



#### Messy Solution

Requires re-entry into sterile field  
Loose fibers



#### Closed System

Sealed to maintain closed system

## Foley Catheter Trays

2000ml drain bag, Sterile, 10 each per case

### Silvertouch 100% Silicone Foley Catheter (Latex-Free)

Item Number	Description
DYND140214	14Fr 5cc, Drain bag with Anti-Reflux Device
DYND140216	16Fr 5cc, Drain bag with Anti-Reflux Device
DYND140218	18Fr 5cc, Drain bag with Anti-Reflux Device

### 100% Silicone Foley Catheter (Latex-Free)

Item Number	Description
DYND11519	16Fr 5cc, Drain bag with Anti-Reflux Device
DYND11520	18Fr 5cc, Drain bag with Anti-Reflux Device

### Lubricious Hydrophilic Latex Foley Catheter

Item Number	Description
DYND11960	16Fr 5cc, Drain bag with Anti-Reflux Device

DYND11903	16Fr 5cc, Drain bag with Anti-Reflux Tower
DYND11908	18Fr 5cc, Drain bag with Anti-Reflux Tower

### Silicone-Elastomer Coated Latex Foley Catheter

Item Number	Description
DYND11855	14Fr 5cc, Drain bag with Anti-Reflux Device
DYND11860	16Fr 5cc, Drain bag with Anti-Reflux Device
DYND11865	18Fr 5cc, Drain bag with Anti-Reflux Device

DYND11003	16Fr 5cc, Drain bag with Anti-Reflux Tower
DYND11008	18Fr 5cc, Drain bag with Anti-Reflux Tower

## Urine Meter Trays

400ml Meter w/2500ml drain bag, Sterile,  
10 each per case

### Silvertouch 100% Silicone Foley Catheter (Latex-Free)

Item Number	Description
DYND140414	14Fr 5cc
DYND140416	16Fr 5cc
DYND140418	18Fr 5cc

### 100% Silicone Foley Catheter (Latex-Free)

Item Number	Description
DYND1516SIL	16Fr 5cc
DYND1518SIL	18Fr 5cc

### Lubricious Hydrophilic Latex Foley Catheter

Item Number	Description
DYND1516LUB	16Fr 5cc
DYND1518LUB	18Fr 5cc

### Silicone-Elastomer Coated Latex Foley Catheter

Item Number	Description
DYND1516SEC	16Fr 5cc
DYND1518SEC	18Fr 5cc



## Add-A-Cath Foley Trays (Latex-Free)

2000ml drain bag, Sterile, 10 each per case

Item Number	Description
DYND11820	10cc Inflation Syringe, Drain Bag with Anti-Reflux Device
DYND11050	30cc Inflation Syringe, Drain Bag with Anti-Reflux Device
DYND11090	10cc Inflation Syringe, Drain Bag with Anti-Reflux Tower
DYND11095	30cc Inflation Syringe, Drain Bag with Anti-Reflux Tower



## Insertion Trays

Medline Foley Catheter Insertion Trays provide the freedom and flexibility to choose the drain bag and foley catheter to fit the needs of your patient. Sterile. 20 each per case.

All trays are **latexfree** PRODUCT

	DYNC 1810	DYNC 1814	DYNC 1815	DYND 10160	DYND 10200	DYND 18100	DYND 10150*
Peel-lid tray	*	*	*				
CSR wrapped tray				*	*	*	*
10cc syringe	*			*		*	*
30cc syringe		*	*		*		
3-pack PVP	*		*	*	*	*	*
3-pack BZK		*					
Aloetouch gloves	*	*	*	*	*	*	*
Specimen Container				*			*

All Trays include:  
 Graduated basin, Aloe 3G gloves,  
 underpad, fenestrated drape, sterile  
 water inflation syringe and lube jelly.



\*DYND10150 also contains a 16FR 5cc Silicone-Elastomer Latex Foley Catheter



## Leg Bags

- Pre-attached leg straps
- Anti-Reflux Valve
- Sterile fluid pathway
- 48 each per case
- All items below are **latexfree** PRODUCT

Item Number	Description
DYND12574	600ml, twist valve drainage port
DYND12575	DYND12574 with extension tubing
DYND12578	1100ml, twist valve drainage port
DYND12579	DYND12578 with extension tubing
DYND12584	600ml, Fabric-Back, exclusive Slide-Tap drainage port
DYND12588	1100ml, Fabric-Back, exclusive Slide-Tap drainage port
DYND12550	18" Extension Tube, 50 ea/cs
DYND12200	Catheter Plug/Protector Cap, 100 ea/cs

## Drain Bags

- Reinforced hanger with built-in anti-kink drain tube guide
- 50" of 11/32" STAR tubing facilitates drainage of thick, fibrous urine while resisting kinking
- Needleless sample port (accessible with either a blunt cannula or needle)
- Adjustable bed sheet clamp secures tubing
- 2000ml capacity
- 20 each per case
- All bags are **latexfree**

## Slide Tap

Item Number	Description
DYND15207	Ultimate Anti-Reflux Device, exclusive Slide-Tap drainage port, infection control insert. Packaged O.R. sterile in Tyvek peel pouch.
DYND15205	Anti-Reflux Tower, exclusive Slide-Tap drainage port. Sterile fluid pathway.
DYNC1674	Anti-Reflux Valve, exclusive Slide-Tap drainage port. Sterile fluid pathway.

## Metal Clamp

Item Number	Description
DYND15210	Ultimate Anti-Reflux Device, drainage port with metal clamp, infection control insert. Packaged O.R. sterile in Tyvek peel pouch.
DYND15203	Anti-Reflux Tower, drainage port with metal clamp. Sterile fluid pathway.
DYND15405	Continuous Bladder Irrigation Drain Bag. 4000ml capacity, Anti-Reflux Tower, drainage port with metal clamp. Packaged O.R. sterile in Tyvek peel pouch.



## Urine Meter

- 400ml capacity dual chamber meter with 2500ml drain bag
- 60" of 9/32" STAR tubing
- O.R. sterile in Tyvek peel pouch
- 10 each per case
- Urine meter is **latexfree**

Item Number	Description
DYND17403	400ml Meter with 2500ml Drain Bag



latexfree  
PRODUCT

## Silvertouch Foley Catheter

100% Silicone Foley Catheters incorporate the power of silver through a patented process that binds silver ions to the catheter's inner and outer lubricious hydrophilic coating. Sterile. 10 each per case.

Item Number	Size	Item Number	Size
<b>2-way Foley, 5cc</b>		<b>2-way Foley, 30cc</b>	
DYND141014	14Fr	DYND141216	16Fr
DYND141016	16Fr	DYND141218	18Fr
DYND141018	18Fr	DYND141220	20Fr
DYND141020	20Fr		



latexfree  
PRODUCT

## 100% Silicone Foley Catheter

100% Silicone Foley Catheters are completely inert to minimize tissue irritation during extended periods of indwelling use. Pliable catheter with a firm tip eases insertion. Features large, smooth eyes for maximum comfort and drainage and symmetrical balloon for added reliability. Sterile. 10 each per case.

Item Number	Size	Item Number	Size
<b>2-way Foley, 3cc</b>		<b>2-way Foley, 30cc</b>	
DYND11552	6Fr	DYND11532	16Fr
DYND11553	8Fr	DYND11533	18Fr
DYND11554	10Fr	DYND11534	20Fr
<b>2-way Foley, 5cc</b>		DYND11535	22Fr
DYND11500	12Fr	DYND11536	24Fr
DYND11501	14Fr		
DYND11502	16Fr		
DYND11503	18Fr		
DYND11504	20Fr		
DYND11505	22Fr		
DYND11506	24Fr		

**For Pre-Connected Closed System Trays, see page 11.**

## Lubricious Hydrophilic Latex Foley Catheter

The lubricious hydrophilic coating on both the inside and outside of the catheter absorbs body fluids while in place to form a protective, slippery barrier between catheter and urethral tissue, resulting in less irritation. Sterile.

12 each per case.


Item Number	Size	Item Number	Size
<b>2-way Foley, 5cc</b>			
DYND11912	12Fr	DYND11932	12Fr
DYND11914	14Fr	DYND11934	14Fr
DYND11916	16Fr	DYND11936	16Fr
DYND11918	18Fr	DYND11938	18Fr
DYND11920	20Fr	DYND11940	20Fr
DYND11922	22Fr	DYND11942	22Fr
DYND11924	24Fr	DYND11944	24Fr
DYND11926	26Fr	DYND11946	26Fr
DYND11928	28Fr	DYND11948	28Fr
DYND11930	30Fr	DYND11950	30Fr



## Silicone-Elastomer Coated Latex Foley Catheter

Innovative manufacturing process provides the maximum amount of silicone coating available. Silicone-elastomer coating resists encrustation. Sterile. 12 each per case.

Item Number	Size	Item Number	Size
<b>2-way Foley, 3cc</b>			
DYND11706	6Fr	DYND11772	12Fr
DYND11708	8Fr	DYND11774	14Fr
DYND11710	10Fr	DYND11776	16Fr
<b>2-way Foley, 5cc</b>			
DYND11752	12Fr	DYND11778	18Fr
DYND11754	14Fr	DYND11780	20Fr
DYND11756	16Fr	DYND11782	22Fr
DYND11758	18Fr	DYND11784	24Fr
DYND11760	20Fr	DYND11786	26Fr
DYND11762	22Fr	DYND11788	28Fr
DYND11764	24Fr	DYND11790	30Fr
DYND11766	26Fr	<b>3-way Foley, 30cc</b>	
DYND11768	28Fr	DYND11800	16Fr
DYND11770	30Fr	DYND11801	18Fr
<b>2-way Foley, 5cc, Coude Tip</b>			
DYND11214	14Fr	DYND11802	20Fr
DYND11216	16Fr	DYND11803	22Fr
DYND11218	18Fr	DYND11804	24Fr
DYND11220	20Fr	DYND11805	26Fr
DYND11222	22Fr		
DYND11224	24Fr		





## From CDC Guideline for Prevention of Catheter-associated Urinary Tract Infection:

“Indwelling catheters should be properly secured after insertion to prevent movement and urethral traction”.

Elastic comfort strap stabilizes the catheter without restricting patient movement. Features a Velcro locking system that is quick and easy to engage. 10 each per box.

Item Number	Description
DYND16900	2" x 22" fits leg with circumference up to 29"
DYND16905	2" x 51" fits leg with circumference up to 65"

Foam catheter strap holds catheter securely in place. 12 each per box.

Item Number	Description
DYND16800	29" x 1-1/2" fits leg with circumference up to 29"

Adhesive patch features a Velcro strap to secure catheter in place. 50 each per case.

Item Number	Description
DYND16700	2" x 3" fits any patient

## Drain Bag Holder

Made from durable poly/cotton. Attaches to wheelchair frame with 100% polyester straps. Three separate snap positions allow easy adjustment.



Item Number	Description	Pkg
MDT825150	Drain bag holder	each

All items listed on this page are **latexfree** PRODUCT

## Catheter Care Tray


Medline catheter care trays are packaged in a peel-lid, double-compartment tray and include the following:

- 1 pair wrapped Aloetouch 3G gloves
- 1 PVP ointment
- 3 PVP swabsticks
- 1 4x4 gauze pad
- 1 PVP prep pad
- 1 Underpad
- Sterile
- 20 each per case

Item Number	Description
DYND12110	Catheter Care Tray


## Cysto/Irrigation Tubing Sets

- 90" length
- Roller clamp
- Drip chamber
- Sterile
- 20 each per case




Item Number	Description
DYND19120	Cysto/TUR Irrigation Set with .188" I.D. tubing
DYND19122	"Y" Type Cysto/TUR Irrigation Set with .281" I.D. tubing, 2 Roberts Clamps above drip chamber
DYND19124	4 Lead Irrigation Set with .281" I.D. tubing, 4 Roberts Clamps above drip chamber





RCH36302



RCH33302



DYND12301

## Male External Catheters

### Self-Adhering 100% Silicone Latex-Free

- **Wide Band®** allows for maximum security, maximum comfort with nonsensitizing adhesive on nearly the entire sheath wall, helping prevent urine migration and build-up.
- **UltraFlex®** has built-in adhesive band for patient comfort and easy application.

### Adhesive Tape Latex

- **Exo-Cath** is a flexible latex catheter packed with double sided adhesive tape to secure catheter.

Item Number	Description	Pkg
RCH36301	<b>Wide Band</b> , Small (25 mm)	30 ea/cs
RCH36302	<b>Wide Band</b> , Medium (29 mm)	30 ea/cs
RCH36303	<b>Wide Band</b> , Intermediate (32 mm)	30 ea/cs
RCH36304	<b>Wide Band</b> , Large (36 mm)	30 ea/cs
RCH36305	<b>Wide Band</b> , Extra-Large (41 mm)	30 ea/cs
RCH33301	<b>UltraFlex</b> , Small (25 mm)	30 ea/cs
RCH33302	<b>UltraFlex</b> , Medium (29 mm)	30 ea/cs
RCH33303	<b>UltraFlex</b> , Intermediate (32 mm)	30 ea/cs
RCH33304	<b>UltraFlex</b> , Large (36 mm)	30 ea/cs
RCH33305	<b>UltraFlex</b> , Extra-Large (41 mm)	30 ea/cs
DYND12303	<b>Exo-Cath</b> , Small (25 mm)	25 ea/cs
DYND12302	<b>Exo-Cath</b> , Medium (30 mm)	25 ea/cs
DYND12301	<b>Exo-Cath</b> , Large (35 mm)	25 ea/cs

## Infant/Pediatric Products

**Infant Urine Collector.** For male or female infant use. Adhesive is gentle on skin and painlessly removed, yet holds securely and folds to form a spill-proof specimen transport. 5 oz (150ml). Latex-free, 50 each per box.


Item Number	Description
MDS190510	Sterile
MDS190505	Non-sterile

**Speci-Cath Kit.** 9" Catheter pre-inserted into specimen collection tube to ensure sterile technique. Also includes identification label and prep supplies (gloves, PVP swabsticks, lubricating jelly). Latex-free. Sterile. 25 each per case.

Item Number	Description
DYND10820	5Fr
DYND10815	8Fr

**Intermittent Urethral Catheters.** Sterile.

Item Number	Description
DYND10712	6Fr Vinyl, 10", latex-free, 30 ea/cs
DYND10713	8Fr Vinyl, 10", latex-free, 30 ea/cs
DYND10714	10Fr Vinyl, 10", latex-free, 30 ea/cs
DYND13508	8Fr Red Rubber latex, 12 ea/cs
DYND13510	10Fr Red Rubber latex, 12 ea/cs




**2-way Foley Catheter with 3cc Balloon.** Sterile. 100% Silicone, Latex-free, 10 each per case.

Item Number	Description
DYND11552	6Fr
DYND11553	8Fr
DYND11554	10Fr

Silicone-Elastomer Coated Latex, 12 each per case.

Item Number	Description
DYND11706	6Fr
DYND11708	8Fr
DYND11710	10Fr



## Clear Vinyl Catheters

Smooth tip, comfort eyes, funnel end. Sterile. Latex-free. Packaged 30 each per case.

### Male/Universal Self Catheter 16"length


Item Number	Description
DYND10721	10Fr
DYND10722	12Fr
DYND10700	14Fr (100 ea/cs)
DYND10723	14Fr
DYND10724	16Fr
DYND10725	18Fr

### Female Self Catheter 6"length

Item Number	Description
DYND10701	10Fr
DYND10702	12Fr
DYND10703	14Fr

### Pediatric Self Catheter 10"length

Item Number	Description
DYND10712	6Fr
DYND10713	8Fr
DYND10714	10Fr



## Red Rubber Latex Catheters

Suitable for male or female catheterization. Smooth tip, two oval eyes and funnel end. 16" length and radiopaque. Sterile.

### Robinson Style

Item Number	Description	Pkg
DYND13508	8Fr	12 ea/cs
DYND13510	10Fr	12 ea/cs
DYND13512	12Fr	12 ea/cs
DYND13514	14Fr	12 ea/cs
DYND13515	15Fr	50 ea/cs
DYND13516	16Fr	12 ea/cs
DYND13518	18Fr	12 ea/cs
DYND13520	20Fr	12 ea/cs
DYND13522	22Fr	12 ea/cs

### Coude-Tip

Item Number	Description	Pkg
DYND13614	14Fr	10 ea/cs

### Red Poly-Cath Latex-Free

Item Number	Description	Pkg
DYND13714	14Fr	10 ea/cs



## Urethral Catheterization Trays

Medline Urethral Catheterization Trays provide you with all of the supplies required for catheterization. Sterile. 20 each per case.



	DYNC 1820	DYNC 1816	DYND 10300	DYND 10305	DYND 10350	DYND 10500
Peel-lid tray	*	*				
CSR wrapped tray			*	*	*	*
Red Rubber Catheter		15Fr			15Fr	
Vinyl Catheter	14Fr		14Fr			14Fr
PVP Swabsticks 3/Pk	*	*	*	*	*	
BZK Swabsticks 3/Pk						*
Aloetouch 3G gloves	*	*	*	*	*	*

All Trays include:  
 graduated basin, Aloe 3G gloves, underpad, fenestrated drape, swabsticks, lube jelly, & specimen container w/label



## My-Cath™ Touch-Free Self Catheter System

The self-catheterization kit for active and independent patients. Features include:

- **Closed System Design** – Special design reduces mess, eases disposal and reduces contamination risk by allowing a touch-free catheterization.
- **Pre-lubricated** – The catheter and introducer are coated with a water soluble lubricant for ease of insertion.
- **Catheter Guide** – The guide allows for simple, protected control of catheter during insertion.
- **1500ml Collection Bag** – The urine collection bag eliminates the need for a separate container, reducing the chance of spillage.

**My-Cath Kit includes:** Collection Bag (as shown above) Moisture-proof underpad, Wrapped gloves, PVP swabsticks (3/pk). Sterile.

Item Number	Description	Pkg
DYND10440	My-Cath Kit	60 ea/cs



latexfree  
PRODUCT

## Pre-Connected Urethral Catheterization Trays

- 1000ml collection bag
- CSR wrapped, sterile
- 20 each per case



All Trays include:  
Aloe 3G gloves, underpad, fenestrated drape, swabsticks, lube jelly, & specimen container w/label



Item Number	Description
DYND10400	With 14Fr vinyl catheter, graduated basin, latex-free
DYND10402	With 14Fr vinyl catheter, foam tray, latex-free
DYND10407	With 15Fr red rubber latex catheter, foam tray

## Speci-Cath Kits latexfree

Medline's Spec-Cath kits provide an easy way to collect a sterile urine sample. The minimally invasive catheter is pre-inserted into the specimen tube to ensure sterile technique. Available in neonatal, pediatric, female and universal/male sizes. Sterile. 25 each per case.

	DYND 10800	DYND 10805	DYND 10806	DYND 10810	DYND 10815	DYND 10820
Catheter size	8Fr	8Fr	8Fr	8Fr	8Fr	5Fr
Catheter length	6 1/2"	6 1/2"	6 1/2"	14"	9"	9"
15ml specimen tube w/cap	*	*	*	*	*	*
ID label	*	*	*	*	*	*
Aloetouch 3G Gloves		*	*	*	*	*
PVP swabsticks 3/pk		*	*	*	*	*
Lubricating Jelly		*		*	*	*





## Deluxe Urinalysis Container

Urinalysis container with pouring spout, 6 oz./180ml, volume graduated in both oz and ml. Broad base helps minimize tipping. Packaged 500 each per case.

Item Number	Description
DYND30100	Clear Plastic container (no lid).
DYND30102	Clear Plastic container with pre-printed lids for specimen identification.
DYND30103	Polypropylene container with pre-printed lids for specimen identification.

## Specimen Container, Polypropylene

With graduations in both oz and ml. Screw-on lid and patient ID label.

Item Number	Description
DYND30330	4 oz./120ml, packaged individually in poly-bag, sterile fluid pathway, 100 each per case.
DYND30331	4 oz./120ml, packaged individually in peel-pouch, sterile, appropriate for O.R. use, 100 each per case.
DYND30335	4 oz./120ml, packaged clean, nested in sleeve, 25 containers per sleeve, 500 each per case.

All items listed on this page are



## Leak Resistant Specimen Container, Polystyrene

Clear, rigid polystyrene container has a metal lid with continuous-thread and "no-leak" seal. Features graduations in both oz and ml. Pre-attached specimen label acts as a tamper-evident seal.

Item Number	Description
DYND30345	4 oz./120ml, packaged bulk, sterile fluid pathway, 200 each per case.

## Leak Resistant Specimen Container, Polypropylene

For use in pneumatic tubing systems

Continuous threaded lid acceptable for use in pneumatic tubing systems. With graduations in both oz and ml. Screw-on lid and patient ID label with tamper evident seal.

Item Number	Description
DYND30337	3 oz./90ml, packaged individually in peel-pouch, sterile, appropriate for O.R. use, 100 each per case.
DYND30332	3 oz./90ml, packaged bulk, sterile fluid pathway, 400 each per case.
DYND30334	4 oz./120ml, packaged bulk, sterile fluid pathway, 100 each per case.
DYND30339	4 oz./120ml, packaged individually in peel-pouch, sterile, appropriate for O.R. use, 100 each per case.

## Valu-Pak Mid-Stream Collection Kits, Sterile

Medline kits include the following components:

- 4 oz./120ml polypropylene specimen container
- Screw-on lid
- Attached protective collar
- Bilingual instructions
- ID label
- Towelette packet

Item Number	Description
DYND30224	Kit with BZK towelettes, 100 each per case.
DYND30228	Kit with Castile Soap towelettes, 100 each per case.



DYND30224

## Basic Mid-Stream Collection Kits, Sterile

Basic kit contains the following components:

- 4 oz./120ml polypropylene specimen container
- Screw-on lid (poly-bagged)
- Bilingual instructions
- ID label
- Towelette packet

Item Number	Description
DYND30240	Kit with BZK towelettes, 100 each per case.



DYND30240

All items listed on this page are **latexfree** PRODUCT

## Deluxe Mid-Stream Collection Sets, Sterile

Deluxe set contains the following components:

- 4 oz./120ml polypropylene specimen container
- Anatomically correct funnel
- Screw-on lid (poly-bagged)
- Bilingual instructions
- ID label
- Ziplock specimen bag
- Towelette packet

Item Number	Description
DYND30211	Kit with Castile Soap towelettes, 20 each per case.
DYND30215	Kit with BZK towelettes, 20 each per case.
DYND30212	Kit with 3 oz./90ml leak-resistant premium container and BZK towelettes, 40 each per case.



DYND30212



## Pathology Containers with Lids

Containers are made of high grade durable polypropylene and can be frozen without cracking. Secure, snap-lock caps minimize leaking.

Item Number	Description	Pkg
DYND34255	8 oz.	250 ea/cs
DYND34260	16 oz.	100 ea/cs
DYND34265	32 oz.	100 ea/cs
DYND34270	64 oz.	50 ea/cs
DYND34275	165 oz.	25 ea/cs



## Specimen Bags

Specimen bags with zip closure and additional document pouch. Includes Biohazard symbol. 1,000 each per case.

Item Number	Description
DYND30260	6"x6" bag with 6"x6" document pouch.
DYND30261	6"x9" bag with 6"x9" document pouch.
DYND30271	8"x8" bag with 8"x8" document pouch.

All items listed on this page are



## Speci-Catch Sputum/Specimen Collector

Ensures a sterile sample with less waste. A 50ml specimen tube with funnel top and cap eliminates the need for transferring. Includes a 3-pk of BZK towelettes and specimen identification label. Sterile.

Item Number	Description
DYND10850	Speci-Catch, 30 each per case



## Stool Specimen Container

Plastic container comes with lid. 250 each per case.

Item Number	Description
DYND30510	Stool Specimen Container



## Polypropylene Urine /Stool Specimen Collector

Fits any standard sized commode. 100 each per case.

Item Number	Description
DYND36600	Translucent, 900ml
DYND36605	Heavy Duty, White, 1200ml



## Commode Specimen Collection Set

700ml capacity container, with snap-on lid, commode stabilizer 100 each per case.

Item Number	Description
DYND36500	Commode Collector



## Irrigation Syringes

A unique barrier flange with "ears" minimize hand slippage and subsequent contamination of barrel or solution container. The 2oz./60ml syringe barrel features large, raised, easy-to-read graduations which are calibrated in both oz and ml. Note: All packaged sterile except DYNC7061 and DYNC7062.

**Contro-Piston Syringe** with thumb ring which is designed to minimize hand slippage and prevent contamination

	Tray	Pad	Sq. Jar	Drape	Luer Tip Protector	Wrapped	Lidding	Qty/Cs
DYNC2303*	Rigid	Alcohol	*	*		*	Tyvek	20
DYNC7061	Rigid		*		*		Tyvek	30
DYNC7062	Soft		*		*		Paper	30
DYND20300	Rigid	PVP	*	*		*	None	20
DYND20302	Rigid	Alcohol	*	*		*	Tyvek	20
DYND20310	Rigid		*			*	Tyvek	20
DYND20325	Soft				*	*	Paper	50
DYND20327	None					*	Bulk	50

\*100ml 0.9% NaCl

**Contro-Bulb Syringe** has soft, pliable, textured Kraton-like polymer bulb for maximum control of drawn or expelled volume. Easy to grip and use with one hand.

	Tray	Pad	Sq. Jar	Drape	Luer Tip Protector	Wrapped	Lidding	Qty/Cs
DYND20100	Rigid	PVP	*	*		*	None	20
DYND20102	Rigid	Alcohol	*	*		*	Tyvek	20
DYND20110	Rigid		*			*	Tyvek	20
DYND20124	Soft					*	Paper	50
DYND20125	Soft					*	Tyvek	50
DYND20127	None					*	Bulk	50



## Feeding Syringe with IV Pole Bag

Feeding syringe comes complete with luer tip adapter and I.V. Pole bag. Single patient use. Non-sterile.

Item Number	Description	Pkg
DYND70631	Piston syringe	30 ea/cs
DYND20335	Piston syringe and jar	30 ea/cs
DYND70651	Bulb syringe	30 ea/cs

All items listed on this page are

**latexfree**  
PRODUCT

## Bard Competitive Cross Reference

Bard Item	Medline Item
000076	DYND12200
35280	DYND20124
35280	DYND20127
35630	DYND10820
35640	DYND10815
35720	DYND10805
35820	DYND10806
38470	DYND20325
38470	DYND20327
120616	DYND13614
150501	DYND12303
150502	DYND12302
150503	DYND12301
150719	DYND12574
150719	DYND12579
150732	DYND12575
150732	DYND12584
150819	DYND12578
150832	DYND12588
153304	DYND17403
153504	DYND15207
153504	DYND15210
153509	DYND15405
154102	DYND15203
154102	DYND15205
165808	DYND11553
165810	DYND11554
165812	DYND11500
165814	DYND11501
165816	DYND11502
165818	DYND11503
165820	DYND11504
165822	DYND11505
165824	DYND11506
166816	DYND11532
166818	DYND11533
166820	DYND11534
166822	DYND11535
166824	DYND11536
277708	DYND13508
277710	DYND13510
277712	DYND13512
277714	DYND13514
277716	DYND13516
277718	DYND13518
277720	DYND13520
277722	DYND13522

Bard Item	Medline Item
420614	DYND10703
421010	DYND10714
421610	DYND10721
421612	DYND10722
421614	DYND10723
750101	DYND20100
750107	DYND20102
750110	DYND20110
750307	DYND20302
750310	DYND20310
771114	DYND10400
771114	DYND10402
771215	DYND10407
772100	DYNC1820
772100	DYND10300
772100	DYND10305
772415	DYNC1816
772415	DYND10350
772415	DYND18350
772417	DYND18200
782102	DYND10200
792100	DYND10160
842801	DYND30215
842803	DYND30211
842904	DYND30228
882816	DYND1516SEC
897416	DYND11519
897418	DYND11520
899016	DYND11003
899018	DYND11008
899616	DYND11903
899618	DYND11908
899716	DYND11960
900014	DYND140214
900016	DYND140216
900018	DYND140218
902116	DYND1516SIL
902118	DYND1518SIL
902816	DYND1516LUB
902818	DYND1518LUB
903114	DYND140414
903116	DYND140416
903118	DYND140418
907300	DYND11090
0100L14	DYND11214
0100L16	DYND11216
0100L18	DYND11218

Bard Item	Medline Item
0100L20	DYND11220
0165L12	DYND11912
0165L14	DYND11914
0165L16	DYND11916
0165L18	DYND11918
0165L20	DYND11920
0165L22	DYND11922
0165L24	DYND11924
0166L12	DYND11932
0166L14	DYND11936
0166L18	DYND11938
0167L16	DYND11800
0167L18	DYND11801
0167L20	DYND11802
0167L22	DYND11803
0167L24	DYND11804
0167L26	DYND11805
165V12S	DYND11752
165V14S	DYND11754
165V16S	DYND11756
165V18S	DYND11758
165V20S	DYND11760
165V22S	DYND11762
165V24S	DYND11764
165V26S	DYND11766
165V28S	DYND11768
165V30S	DYND11770
166V12S	DYND11772
166V14S	DYND11774
166V16S	DYND11776
166V18S	DYND11778
166V20S	DYND11780
166V22S	DYND11782
166V24S	DYND11784
166V26S	DYND11786
166V28S	DYND11788
166V30S	DYND11790
1758SI14	DYND141014
1758SI16	DYND141016
1758SI18	DYND141018
1758SI20	DYND141020
1768SI16	DYND141216
1768SI18	DYND141218
1768SI20	DYND141220
4A4194	DYND12550
4A5144	DYND10440



## Kendall Competitive Cross Reference

Kendall Item	Medline Item
1416	DYND11932
1418	DYND11936
1612	DYND11805
1614	DYND11912
1616	DYND11914
1618	DYND11916
1620	DYND11918
1622	DYND11920
1624	DYND11922
2006	DYND1518SEC
2008	DYND1516SIL
2018	DYND1516LUB
2552	DYND10702
3305	DYND10300
3411	DYND10805
3433	DYND12574
6014	DYND11903
6080	DYND11908
6140	DYND11095
6210	DYND11519
6255	DYND11020
6300	DYND15205
7000	DYND1518SIL
8206	DYND15203
8256	DYND11050
8946	DYND11520
8948	DYND11003
22000	DYND30212
24000	DYND30228
25000	DYND30215
25018	DYND30224
67000	DYND20124
67800	DYND20100
67803	DYND70631
68000	DYND20310
68800	DYND20300
68803	DYNC7062
75000	DYND10350
75030	DYNC1816
75035	DYND10305
145514	DYND15405
145516	DYND12575
400608	DYND10714
400610	DYND10721
400614	DYND10722
400616	DYND10723
400618	DYND10724

Kendall Item	Medline Item
402712	DYND11710
402714	DYND11752
402716	DYND11754
402718	DYND11756
402720	DYND11758
402722	DYND11760
402724	DYND11762
402726	DYND11764
402728	DYND11766
402730	DYND11768
403712	DYND11770
403714	DYND11772
403716	DYND11774
403718	DYND11776
403720	DYND11778
403722	DYND11780
403724	DYND11782
403726	DYND11784
403728	DYND11786
403730	DYND11788
8884732000	DYND12301
8887600032	DYND10402
8887600776	DYND18100
8887601220	DYND10200
8887603085	DYND11552
8887603101	DYND11553
8887605122	DYND11938
8887605148	DYND11500
8887605163	DYND11501
8887605189	DYND11502
8887605205	DYND11503
8887605221	DYND11504
8887605247	DYND11505
8887630187	DYND11532
8887630203	DYND11533
8887630229	DYND11534
8887630245	DYND11535
8887654088	DYND11706
8887654104	DYND11708
8887660085	DYND10725
8887660101	DYND13508
8887660127	DYND13510
8887660143	DYND13512
8887660168	DYND13515
8887660184	DYND13516
8887660218	DYND13518
8887671660	DYND11855

Kendall Item	Medline Item
8887671686	DYND11860
8887689167	DYND11220
8887689183	DYND11800
8887689209	DYND11801
8887689225	DYND11802
8887689241	DYND11803
8887689266	DYND11804
1614C	DYND11790
1616C	DYND11214
1618C	DYND11216
1620C	DYND11218
314V2	DYND10407
3450LF	DYND10400
605148IC	DYND140418
605163IC	DYND141014
605189IC	DYND141016
605205IC	DYND141018
6146IC	DYND140214
6148IC	DYND140216
630187IC	DYND141216
630203IC	DYND141218
7004IC	DYND140218
7006IC	DYND140414
7008IC	DYND140416

Item	Page	HCPSC	Latex-Free
DYNC1674	13	A4357	Latex-Free
DYNC1810	12	A4310	Latex-Free
DYNC1814	12	A4310	Latex-Free
DYNC1815	12	A4310	Latex-Free
DYNC1816	18	A4353	
DYNC1820	18	A4353	Latex-Free
DYNC2303	23	A4320	Latex-Free
DYNC7061	23	A4320	Latex-Free
DYNC7062	23	A4320	Latex-Free
DYND10150	12	A4311	
DYND10160	12	A4310	Latex-Free
DYND10200	12	A4310	Latex-Free
DYND10300	18	A4353	Latex-Free
DYND10305	18	A4310	Latex-Free
DYND10350	18	A4353	
DYND10400	19	A4353	Latex-Free
DYND10402	19	A4353	Latex-Free
DYND10407	19	A4353	
DYND10440	19	A4353	Latex-Free
DYND10500	18	A4353	Latex-Free
DYND10700	18	A4351	Latex-Free
DYND10701	18	A4351	Latex-Free
DYND10702	18	A4351	Latex-Free
DYND10703	18	A4351	Latex-Free
DYND10712	17,18	A4351	Latex-Free
DYND10713	17,18	A4351	Latex-Free
DYND10714	17,18	A4351	Latex-Free
DYND10721	18	A4351	Latex-Free
DYND10722	18	A4351	Latex-Free
DYND10723	18	A4351	Latex-Free
DYND10724	18	A4351	Latex-Free
DYND10725	18	A4351	Latex-Free
DYND10800	19	A9270	Latex-Free
DYND10805	19	A9270	Latex-Free
DYND10806	19	A9270	Latex-Free
DYND10810	19	A9270	Latex-Free
DYND10815	17,19	A9270	Latex-Free
DYND10820	17,19	A9270	Latex-Free
DYND10850	22		Latex-Free
DYND11003	11	A4314	
DYND11008	11	A4314	
DYND11050	11	A4354	Latex-Free
DYND11090	11	A4354	Latex-Free
DYND11095	11	A4354	Latex-Free
DYND11214	15	A4340	
DYND11216	15	A4340	
DYND11218	15	A4340	
DYND11220	15	A4340	
DYND11222	15	A4340	
DYND11224	15	A4340	
DYND11500	14	A4344	Latex-Free
DYND11501	14	A4344	Latex-Free
DYND11502	14	A4344	Latex-Free
DYND11503	14	A4344	Latex-Free
DYND11504	14	A4344	Latex-Free
DYND11505	14	A4344	Latex-Free
DYND11506	14	A4344	Latex-Free
DYND11519	11	A4315	Latex-Free
DYND11520	11	A4315	Latex-Free

Item	Page	HCPSC	Latex-Free
DYND11532	14	A4344	Latex-Free
DYND11533	14	A4344	Latex-Free
DYND11534	14	A4344	Latex-Free
DYND11535	14	A4344	Latex-Free
DYND11536	14	A4344	Latex-Free
DYND11552	14,17	A4344	Latex-Free
DYND11553	14,17	A4344	Latex-Free
DYND11554	14,17	A4344	Latex-Free
DYND11706	15,17	A4338	
DYND11708	15,17	A4338	
DYND11710	15,17	A4338	
DYND11752	15	A4338	
DYND11754	15	A4338	
DYND11756	15	A4338	
DYND11758	15	A4338	
DYND11760	15	A4338	
DYND11762	15	A4338	
DYND11764	15	A4338	
DYND11766	15	A4338	
DYND11768	15	A4338	
DYND11770	15	A4338	
DYND11772	15	A4338	
DYND11774	15	A4338	
DYND11776	15	A4338	
DYND11778	15	A4338	
DYND11780	15	A4338	
DYND11782	15	A4338	
DYND11784	15	A4338	
DYND11786	15	A4338	
DYND11788	15	A4338	
DYND11790	15	A4338	
DYND11800	15	A4346	
DYND11801	15	A4346	
DYND11802	15	A4346	
DYND11803	15	A4346	
DYND11804	15	A4346	
DYND11805	15	A4346	
DYND11820	11	A4354	Latex-Free
DYND11855	11	A4314	
DYND11860	11	A4314	
DYND11865	11	A4314	
DYND11903	11	A4314	
DYND11908	11	A4310	
DYND11912	15	A4338	
DYND11914	15	A4338	
DYND11916	15	A4338	
DYND11918	15	A4338	
DYND11920	15	A4338	
DYND11922	15	A4338	
DYND11924	15	A4338	
DYND11926	15	A4338	
DYND11928	15	A4338	
DYND11930	15	A4338	
DYND11932	15	A4338	
DYND11934	15	A4338	
DYND11936	15	A4338	
DYND11938	15	A4338	
DYND11940	15	A4338	
DYND11942	15	A4338	

Item	Page	HCPCS	Latex-Free
DYND11944	15	A4338	
DYND11946	15	A4338	
DYND11948	15	A4338	
DYND11950	15	A4338	
DYND11960	11	A4314	
DYND12110	16	A4310	Latex-Free
DYND12200	12		Latex-Free
DYND12301	17	A4349	
DYND12302	17	A4349	
DYND12303	17	A4349	
DYND12550	12	A4331	Latex-Free
DYND12574	12	A4358	Latex-Free
DYND12575	12		Latex-Free
DYND12578	12	A4358	Latex-Free
DYND12579	12		Latex-Free
DYND12584	12	A4358	Latex-Free
DYND12588	12	A4358	Latex-Free
DYND13508	17,18	A4351	
DYND13510	17,18	A4351	
DYND13512	18	A4351	
DYND13514	18	A4351	
DYND13515	18	A4351	
DYND13516	18	A4351	
DYND13518	18	A4351	
DYND13520	18	A4351	
DYND13522	18	A4351	
DYND13614	18	A4352	
DYND13714	18	A4351	Latex-Free
DYND140214	9,11		Latex-Free
DYND140216	9,11		Latex-Free
DYND140218	9,11		Latex-Free
DYND140414	9,11		Latex-Free
DYND140416	9,11		Latex-Free
DYND140418	9,11		Latex-Free
DYND141014	9,14		Latex-Free
DYND141016	9,14		Latex-Free
DYND141018	9,14		Latex-Free
DYND141020	9,14		Latex-Free
DYND141216	9,14		Latex-Free
DYND141218	9,14		Latex-Free
DYND141220	9,14		Latex-Free
DYND1516LUB	11	A4314	
DYND1516SEC	11	A4314	
DYND1516SIL	11	A4315	Latex-Free
DYND1518LUB	11	A4314	
DYND1518SEC	11	A4314	
DYND1518SIL	11	A4315	Latex-Free
DYND15203	13	A4357	Latex-Free
DYND15205	13	A4357	Latex-Free
DYND15207	13	A4357	Latex-Free
DYND15210	13	A4357	Latex-Free
DYND15405	13	A4357	Latex-Free
DYND16700	16	A4333	Latex-Free
DYND16800	16	A4334	Latex-Free
DYND16900	16	A4334	Latex-Free
DYND16905	16	A4334	Latex-Free
DYND17403	13	A4357	Latex-Free
DYND18100	12	A4310	Latex-Free
DYND19120	16		Latex-Free

Item	Page	HCPCS	Latex-Free
DYND19122	16		Latex-Free
DYND19124	16		Latex-Free
DYND20100	23	A4320	Latex-Free
DYND20102	23	A4320	Latex-Free
DYND20110	23	A4320	Latex-Free
DYND20124	23	A4322	Latex-Free
DYND20125	23	A4322	Latex-Free
DYND20127	23	A4322	Latex-Free
DYND20300	23	A4320	Latex-Free
DYND20302	23	A4320	Latex-Free
DYND20310	23	A4320	Latex-Free
DYND20325	23	A4322	Latex-Free
DYND20327	23	A4322	Latex-Free
DYND20335	23		Latex-Free
DYND30100	20		Latex-Free
DYND30102	20		Latex-Free
DYND30103	20		Latex-Free
DYND30211	21		Latex-Free
DYND30212	21		Latex-Free
DYND30215	21		Latex-Free
DYND30224	21		Latex-Free
DYND30228	21		Latex-Free
DYND30240	21		Latex-Free
DYND30260	22		Latex-Free
DYND30261	22		Latex-Free
DYND30271	22		Latex-Free
DYND30330	20		Latex-Free
DYND30331	20		Latex-Free
DYND30332	20		Latex-Free
DYND30334	20		Latex-Free
DYND30335	20		Latex-Free
DYND30337	20		Latex-Free
DYND30339	20		Latex-Free
DYND30345	20		Latex-Free
DYND30510	22		Latex-Free
DYND34255	22		Latex-Free
DYND34260	22		Latex-Free
DYND34265	22		Latex-Free
DYND34270	22		Latex-Free
DYND34275	22		Latex-Free
DYND36500	22		Latex-Free
DYND36600	22		Latex-Free
DYND70631	23	B4034/5/6	Latex-Free
DYND70651	23		Latex-Free
MDS190505	17	A4328	Latex-Free
MDS190510	17	A4328	Latex-Free
MDT825150	16		Latex-Free
RCH33301	17	A4349	Latex-Free
RCH33302	17	A4349	Latex-Free
RCH33303	17	A4349	Latex-Free
RCH33304	17	A4349	Latex-Free
RCH33305	17	A4349	Latex-Free
RCH36301	17	A4349	Latex-Free
RCH36302	17	A4349	Latex-Free
RCH36303	17	A4349	Latex-Free
RCH36304	17	A4349	Latex-Free
RCH36305	17	A4349	Latex-Free

# What's NEW in Medline Urologicals?



Silvertouch. The Standard for Silver.

See page 8-9 for more information.



**1-800-MEDLINE** | [www.medline.com](http://www.medline.com)

©2006 Medline Industries, Inc. Medline is a registered trademark and Silvertouch is a trademark of Medline Industries, Inc. WideBand and Ultra Flex are registered trademarks of Rochester Medical Corp. MKT206235/LIT184/10M/ZPN6