



MED-CHECKS

STERILIZATION MONITORING PRODUCTS

WHEN ASSURING STERILITY IS YOUR JOB YOU CAN COUNT ON MEDLINE'S QUALITY LINE OF STERILIZATION MONITORING PRODUCTS





Med-Checks Bowie Dick Advanced Steam Sterilizer Tests

The Med-Checks Bowie Dick advanced test offers accurate air removal results without the cumbersome packing material of other tests.

- Detects air leaks in a pre vacuum steam sterilizer at 270°F (132°C) at a minimum 3.5-minute exposure
- Unique test card design does not require packing material for accurate readings
- Reusable holder saves space on your shelf and reduces waste
- Distinctive color change (from purple to green) that is easy to read and not reversible
- Load record log on the back of each test card provides room for all necessary information
- Lead-free ink and free of any toxic heavy metals

Starter Set

MDS200450S

31 Test Cards, 1 Holder

Test cards-refills

MDS200450

30/pk, 120/cs



MDS200450S



MDS200450

Med-Checks Bowie Dick Test Pack, Steam

- Detects air leaks in a pre vacuum steam sterilizer at 270°F (132°C) at a minimum 3.5-minute exposure
- Distinctive color change (from blue to black) that is not reversible
- Lead-free ink and free of any toxic heavy metals

Steam Test Pack

MDS200400

20/cs

Steam Test Sheets

MDS200700

8 1/2" x 11" 50/pk, 600/cs



MDS200400

Sterilization Indicators and Integrators

Med-Checks Steam Indicator Strips

- Distinctive color change (from blue to black) that is not reversible
- For use in gravity and vacuum cycles

MDS200100	4"	250/bx
MDS200100C	4"	250/bx, 4500/cs
MDS200150	8"	250/bx
MDS200150C	8"	250/bx, 2250/cs



MDS200100

Med-Checks Gas Integrator Strips

- Measures the time, temperature, humidity, and EO gas concentration
- Distinctive color change (from yellow to brown) that is not reversible
- Plastic laminated

MDS200200	4" Strips	250/bx
-----------	-----------	--------



MDS200600

Med-Checks Steam Integrator Strips

- Measures the time, temperature and presence of steam in steam sterilizers
- Distinctive color change (from purple to green) that is not reversible
- Plastic laminated

MDS200600	4"	250/bx
MDS200600C	4"	250/bx, 2500/cs

Med-Checks Plasma Indicator Strips

- Confirms the strip has been exposed to hydrogen peroxide during processing
- Distinctive color change (from blue to pink) that is not reversible
- Color change occurs within 30 minutes of exposure to hydrogen peroxide during the diffusion phase of processing

MDS200300	4"	250/bx
-----------	----	--------



MDS200836 and MDS200848

Med-Checks Steam Autoclave Tape

- Distinctive color change (from beige to black) that is not reversible

MDS200836	1" x 60 yds.	36 rolls/cs
MDS200848	3/4" x 60 yds.	48 rolls/cs



Serilization Record Keeping

Med-Checks Sterilization Record Envelope, Steam/Gas

- Preprinted 9" x 12" envelope with storage capacity for Bowie Dick test sheets, sterilization run records and other pertinent data

MDS200900 100/pk

Med-Checks Sterilization Record Card, Steam/Gas

- Test for both steam and EO gas sterilization
- Information, such as sterilization date or type of load can be written on the back of the card
- Distinctive color change (from blue to black) during steam sterilization and (from yellow to brown) during EO sterilization that is not reversible

MDS200800 3 1/2" x 5" 500/bx



MDS200800 and MDS200900

Biological Indicators and Test Packs

Med-Checks Biological Incubator

- 14 vial capacity for use with steam (57°C) or gas (37°C) indicators - Includes Thermometer

MDS200580 Incubator each

Med-Checks Biological Test Pack & Indicator, Steam

- Self-contained biological indicator that can be incubated in any small incubator or in a 3M® incubator
- Distinctive color change (from purple to yellow) will occur if the spores survive

MDS200500 Test Pack 30/cs
(24 hour incubation)

MDS200525 Indicators 100/bx
(24 hour incubation)

Med-Checks Biological Test Pack, Instant Readout, Steam

- Contains an EZ test self-contained biological indicator and an Instant Readout Integrator card
- After sterilization, the Instant Readout Integrator card predicts the outcome of the biological test for immediate load release
- Distinctive color change (from purple to green) Instant Readout Integrator color change
- Biological Indicator will change (from purple to yellow) if the spores survive

MDS200500IR Test Pack 30/cs

Ask about Medline's full-line of Central Sterile products including: sterilization pouches & monitoring products, instrument care & handling, and sterilization containers & accessories.

