

## **Ultrasorbs**™ Premium Disposable Underpad

The most advanced underpad for optimal skin dryness Keeps the bedding and patient dry

**Available** 

**Ultrasorbs Breathable** 

For use on low air loss beds

"We are very impressed with Ultrasorb's performance. It is the only undernad we have seen that wicke away fluid and keens the nation! "We are very impressed with Ultrasorb's performance. It is the only and keeps the patient's wicks away fluid and keeps the patient's at night as underpad we have seen that wicks away for start onen airing at night as underpad we have seen that wicks away fluid and are represented in the patient's at night as underpad we have seen that wicks allowed us to start onen airing at night. underpad we have seen that wicks away fluid and keeps the patient, of our modern of ou skin dry. Ultrasorbs has allowed us to start open-airing at night as

Skin dry. Ultrasorbs has allowed us to start open-airing at night as

part of our updated pressure ulcer prevention program.

Part of our updated pressure ulcer prevention do a change easier and faster

part of our updated pressure ulcer prevention program. part of our updated pressure ulcer prevention program. Our Patients are more comfortable, and the staff can do a change easier and faster are more comfortable, and the staff can do a change easier and the staff can do a change easier and faster are more comfortable, and the staff can do a change easier and faster with briefs." ramme Obura-Wikes, KN

Director of staff development Wound care coordinator, Pauline Obura Wilkes, RN than with briefs."

"Ultrasorbs has saved us over 20% in product cost alone because of We went from using an it's avtraordinary absorbancy and drynass. We went from using an "Ultrasorbs has saved us over 20% in product cost alone because of We went from using an it's extraordinary absorbency and dryness. We went to were also it's extraordinary absorbency are iner one Illrasorbe We went also also average of 3 – 4 undernade to iner one Illrasorbe was average was ave it's extraordinary absorbency and dryness. We went from using an average of 3 – 4 underpads to just one consistent average of 3 the strength and the consistent average of 3 impressed with the strength and the consistent average. average of 3 - 4 underpads to just one Oltrasorns. We wanted average of 3 - 4 underpads to just one consistent quality."

impressed with the strength and the consistent quality."

DON, Skilled Nursing Facility, Sacramento, California

"We standardized all our underpads to Ultrasorbs we standardized all our underpads to Ultrasorbs We have because of it's high quality and absorbency.

because of it's high quality are product weare hereing the product weareness to ultrasorbs. because or it's high quality and absorbency. We have reduced our staff time as well as product usage because reduced on staff time as well as product The net result reduced our start time as well as product usage because

The net result as often. The net result as often.

We are not having to change as often.

We are not having to change as often.

is savings and better care." Materials Manager, Hospital Houston, Texas

As noted in WOCN 2003 Guidelines,

Managing skin moisture is a critical factor in maintaining skin integrity and reducing the incidence of pressure ulcers.

**Ultrasorbs wicks moisture away from** the patient's skin and is super absorbent to provide optimal skin care.



### **Ultrasorbs™ Premium Disposable Underpad**

- Advanced Supercore™ thermo-bonded core delivers superior skin dryness and resident comfort.
- Wicks moisture away from the skin and locks fluid away to provide better skin care and resident dignity.
- Super absorbent | Ultrasorbs absorbs what used to require 3 or more underpads.
- Super strong to reduce tearing and resulting linen changes.
- Superior odor control.

### **Unique Construction**

- Folded over edge for better containment.
- Soft non-woven top sheet across the pad protects skin from contact with the poly backing.
- Ultra absorbent polymer is spread throughout the entire core for maximum absorbency.
- Top sheet is fused to the core so it does not stick to the skin or bunch up while in use.



Advanced SuperCore™ (Thermo Bonded Core)

# Now available in two styles

### Ultrasorbs™

• Waterproof polyethylene backing protects underlying surfaces.

Item No: MSC282100

Size:  $23'' \times 36''$  overall,

21" × 32" absorbent core

Packaging: 70 per case (bags of 10)

### **Ultrasorbs<sup>™</sup> Breathable**

- For use on low air loss beds.
- Spun-bond breathable backing allows air to flow up through the pad while still protecting the bed.

Item No: MSC282500 Size: 23" × 36" overall.

 $21'' \times 32''$  absorbent core

Packaging: 60 per case (bags of 10)

